



Hillrom is a part of Baxter

Baxter



Hillrom™

Welch Allyn® Connex®
Vital Signs Monitor

THE MEDICAL SURGICAL FLOOR YOU'RE USED TO IS CHANGING

Patient deterioration on the MedSurg floor needs a vital signs solution that can monitor continuously, even when you can't. The Welch Allyn® Connex® Vital Signs Monitor can help.



THE RIGHT LEVEL OF MONITORING At your fingertips

Capture standard vitals — like temperature, blood pressure, pulse rate and SpO₂ — in a single, automated solution. Customise your monitoring needs with a streamlined transition from spot check to continuous monitoring.



MONITOR VENTILATION With capnography

Unanticipated respiratory decline can be a challenge. We can help you ensure patient safety with optional EtCO₂ for patients who are at risk of respiratory failure. The Integrated Pulmonary Index (IPI) provides a simplified respiratory status value by combining EtCO₂, respiration rate, pulse rate and SpO₂ information.



TRACK CARDIAC ACTIVITY Using 3- or 5-lead ECG monitoring

Scale up with ECG monitoring to help identify lethal arrhythmia events. Streamline your workflow by adding cardiac assessment to your monitoring routine.



IDENTIFY PATIENTS AT RISK Through early warning scoring protocols

Customise early warning scores to meet your facility's protocols.



SIMPLE

Take a full set¹ of patient vitals and see them on device screen in less than one minute, and electronically send the information to your patient's chart. Operate on an easy-to-use interface for both spot check and continuous workflows. You can switch between monitoring workflows with a streamlined process to meet your patient's monitoring needs.

Capnography, ECG monitoring and continuous SpO₂ are available for patients who require additional attention. Stay focused on patient care with device alarms that alert you to important vital sign changes.

SECURE

Connex Vital Signs solutions are U.S. Department of Defense (DoD) ATO certified, so you can protect patient information on your device and over your network. Securely manage data with encryption, single sign-on, screen lock and server authentication.

CONNECTED

Significantly reduce vital signs documentation errors and give clinicians time back in their day with automated EMR connectivity. Securely send data directly from the bedside to a patient's chart in your EMR so all authorised caregivers can access the most up-to-date information. And with a connection to the Welch Allyn Connex Central Station, clinicians can add a layer of visibility by monitoring patients from a remote, centralised location.²

PRODUCT ORDERING INFORMATION³

Part Number	Description
67MXTX-6	SureBP [®] NIBP, SureTemp [®] Plus Thermometry, Masimo [®] SpO ₂ (Add WiFi: 68MXTX-6)
67NXTX-6	SureBP NIBP, SureTemp Plus Thermometry, Nellcor [®] SpO ₂ (Add WiFi: 68NXTX-6)
68NCTX-6	Capnography (Covidien [®]), SureBP NIBP, SureTemp Plus Thermometry, Nellcor SpO ₂ , Wireless Radio

ACCESSORIES⁴

Part Number	Description
4800-60	Accessory Cable Management mobile stand with basket
4900-60	Accessory Cable Management Stand with basket, extended housing
008-0834-01	GCX wall mount
4900-62	GCX wall mount, extended housing
06000-900	Braun PRO 6000 dock with thermometer
6000-916HS	Barcode scanner with coiled USB
6000-ECG3A	Three-lead ECG module
6000-ECG5A	Five-lead ECG module

SMARTCARE™ SERVICE PROGRAM OPTIONS

Program	Part Number	Program Terms
Protection	S1-CVSM-PRO-PS	3 Years
Protection+	S9-CVSM-PROPLUS-PS	3 Years



Hillrom[™]

To learn more about upgrading your Vitals devices, contact your Hillrom representative or call Hillrom Customer Care on 1800 445 576.

Always read the labels and instructions for use. Contact your Hillrom Representative for the latest versions.

¹ Includes non-invasive blood pressure, SpO₂, heart rate and oral temperature.

² Remote monitoring on ECG is not supported.

³ Sample product ordering information. Some of these configurations are for the US region only, so please contact your Hillrom representative to find the best solution for your facility and monitoring needs.

⁴ Additional accessories are available. Please contact your Hillrom representative to learn more.

Covidien is a registered trademark of Covidien AG. Masimo is a registered trademark of Masimo Corp. Nellcor is a registered trademark of Covidien LP.

Connex Vital Signs solutions are U.S. Department of Defence (DoD) ATO certified*. (*DoD ATO Certification data on Hillrom file.)