



Hillrom™

Welch Allyn FlexiPort® Disposable Blood Pressure Soft Cuffs



Single Patient use blood pressure cuffs can help enhance patient flow and reduce the risk of cross-contamination.

Welch Allyn offers Soft Cuff and EcoCuff® ranges of Disposable Blood Pressure Cuffs to minimise patient to patient transmission of HAI's in clinical settings.



FLEXIPOINT DISPOSABLE SOFT CUFF



SOFT CUFF	ECOCUFF
Large Adult 30cm – 43cm	Large Adult 33cm – 45cm
Adult 25cm – 34cm	Adult 27cm – 38cm
Small Adult 20cm – 26cm	Small Adult 20cm – 28cm
Child 15cm – 21cm	Child 15cm – 21cm



FLEXIPOINT ECOCUFF SINGLE PATIENT USE CUFF



FEATURE	FLEXIPOINT SOFT CUFF	FLEXIPOINT ECOCUFF
Durability	Rated to 500 Inflations	Rated to 100 Inflations Bio-degradable/Recyclable
Appearance	Printed Part Numbers Inverse Colour Coded by Size	Appropriate Sizing Slot
Fastening	Velcro Fastening	
Presentation	Box of 20	Case of a 100

To learn more or call Hillrom Customer Care on 1800 445 576

Always read the labels and instructions for use. Contact your Hillrom Representative for the latest version.

AUTHORISED AUSTRALIAN SPONSOR Welch Allyn Australia Pty Limited
Suite 4.01, 2 Lyonpark Road, Macquarie Park NSW Australia 2113

© 2021 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR120301-AU_r1 02/2021

hillrom.com.au

Welch Allyn FlexiPort® Disposable Blood Pressure Soft Cuffs

FLEXIPOINT DISPOSABLE SOFT CUFF



CODE	DESCRIPTION
SOFT-06	Small Infant, 7-10 cm
SOFT-07	Infant, 9-13 cm
SOFT-08	Small Child, 12-16 cm
SOFT-09	Child, 15-21 cm
SOFT-10	Small Adult, 20-26 cm
SOFT-11	Adult, 25-34 cm
SOFT-11L	Adult Long, 25-34 cm
SOFT-12	Large Adult, 32-43 cm
SOFT-12L	Large Adult Long, 32-43 cm
SOFT-13	Thigh, 40-55 cm

Welch Allyn FlexiPort blood pressure cuffs can help enhance patient flow and reduce the risk of cross-contamination. FlexiPort makes it easier to find the proper-sized cuff for your patient and help prevent miscuffing.



FLEXIPOINT ECOCUFF SINGLE PATIENT USE CUFF



CODE	DESCRIPTION
ECOCUFF-09	Child, 15-21 cm
ECOCUFF-10	Small Adult, 20-28 cm
ECOCUFF-11	Adult, 27-38 cm
ECOCUFF-12	Large Adult, 33-45 cm

One true solution for reducing costs, cross-contamination and environmental impact.



Hillrom™