

Baxter

Hillrom Centrella

SMART+ BED



**THERAPEUTIC
SURFACES**

INNOVATIVE SUPPORT SURFACES HELP MANAGE PRESSURE, SHEAR AND MOISTURE TO DELIVER OPTIMAL WOUND PREVENTION AND HEALING



CENTRELLA SMART+ BED CORE SURFACE

NON-POWERED PRESSURE
REDISTRIBUTION DESIGNED TO
ENHANCE PATIENT COMFORT

- Sloped heel section can help offload pressure in the vulnerable heel area
- Four-way stretch cover helps with patient immersion and envelopment
- RF welded seams help prevent fluid ingress



Multi-Layer Foam Construction

Up to 11% less peak sacral pressure.¹



CENTRELLA SMART+ BED PRO SURFACE

TARGETED AND CUSTOMISED NON-
POWERED PRESSURE REDISTRIBUTION
FOR OPTIMAL PATIENT SUPPORT AND
COMFORT

- Sloped heel section can help offload pressure in the vulnerable heel area
- Viscoelastic foam contours to patient's unique body shape, increasing envelopment, helping to improve patient comfort
- 91 cm – 101 cm Width options support broader range of patients with higher BMIs



Targeted Pressure Redistribution

Up to 21% less peak sacral pressure.¹



Best-In-Class Comfort

Rated 5X more comfortable.²

Proven solutions to help reduce pressure injuries and their associated costs.



IPUP™ SURVEY



HILLROM AND WOCN®
SURFACE ALGORITHM



HILLROM SURFACE
TECHNOLOGIES

= **33%**



CENTRELLA SMART+ BED PRO+ SURFACE

INTEGRATED ADVANCED **MICROCLIMATE** THERAPY HELPS MANAGE HEAT AND MOISTURE IN SECONDS

- Sloped heel section can help offload pressure in the vulnerable heel area
- Provide **microclimate** management on demand with the integrated pump
- Integrated blower design can help improve caregiver efficiency and therapy delivery
- 91 cm – 101 cm Width options support broader range of patients with higher BMIs



Focused Advanced Microclimate Technology

Removes up to 2.5 X more heat and 4 X more moisture.³



Targeted and Customised Pressure Redistribution

Up to 8% less peak sacral pressure.³



CENTRELLA SMART+ BED MAX SURFACE

ADVANCED TECHNOLOGY DELIVERING SUPERIOR PRESSURE REDISTRIBUTION AND **MICROCLIMATE** MANAGEMENT

- Side-sitting detection to help provide stability and support
- Advanced bottoming out protection for patient safety
- Patient comfort adjust for personalised settings



Focused Advanced Microclimate Technology

Removes up to 4 X more heat and 7 X more moisture.²



Weight-Based Pressure Redistribution

Up to 13% less peak sacral pressure.¹



Enhanced Turn Assist

More effective turn in 77% less time.¹

For more information, please contact your Baxter representative
or call Baxter Customer Care on 1800 445 576.

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE

CONTACT YOUR HILLROM REPRESENTATIVE FOR THE LATEST VERSIONS

References

1 Lachenbruch C. 2017. Data on file.

2 Smith K. 2017. Data on file.

3 Hillrom Internal Test Results. Data on file.

4 Hillrom Impact Tracker Case Study 193131. Acute Care Hospital Reduces Pressure Ulcer Rates 33%. Sept. 2015.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

AUTHORISED AUSTRALIAN SPONSOR Hill-Rom Pty. Ltd. 1 Baxter Drive, Old Toongabbie NSW 2146.

hillrom.com.au

Baxter

AU_CS3-230015_r1 07/2023