

Connex[®] Vital Signs Monitor

with Braun ThermoScan[®] PRO 6000 Ear Thermometry Option

When you purchase the Connex Vital Signs Monitor (VSM), you can choose SureTemp[®] Plus oral, pediatric axillary, adult axillary and rectal thermometry PLUS Braun ThermoScan PRO 6000 ear thermometry. You can opt to include one or both thermometry options and both can be fully integrated with the Connex VSM, allowing for direct communication of readings to the Connex VSM screen and the transfer of data to your electronic medical record (EMR).



Plus, the Braun Ear Thermometer can be configured with your Connex VSM so it works best for your facility workflow:



Mount the Braun ThermoScan PRO 6000 directly onto the Welch Allyn accessory management stand.

06000-900
(snaps into stand shelf)



Mount the Braun ThermoScan PRO 6000 directly onto the Welch Allyn economy mobile stand.

06000-900 and 86200
hardware kit



Utilize the tabletop option for the Braun ThermoScan PRO 6000 if you plan to keep the monitor in a stationary location such as triage area, clinic or smaller office.

06000-900
(includes cradle tabletop foot)



Mount the Braun ThermoScan PRO 6000 directly to the wall—mount basket for easy access.

06000-900 and 86100
hardware kit

Braun ThermoScan PRO 6000 Specifications

Performance	Calibration accuracy of ± 0.4 °F (± 0.2 °C) for the thermometer at 95.9 °F - 107.6 °F (35.5 °C - 42.0 °C)
	Response time: 3-7 seconds (Initial start up 3-4 seconds, measurement 2-3 seconds)
Temperature range	Ambient temperature range: 50.0 °F - 104.0 °F (10.0 °C - 40.0 °C)
	Patient temperature range: 68.0 °F - 108.0 °F (20.0 °C - 42.2 °C)
	Display resolution: 0.1 °F or °C (one decimal point displayed)
Dimensions	Weight (without batteries): 3.6 ounces (100 grams)
	Unit dimensions: 5.98" x 1.73" x 1.30" (152 mm x 44 mm x 33 mm)
Power	Rechargeable batteries: Approx. 500 readings per charge and 3 months standby time
Warranty	3-year limited instrument, and 1-year Braun PRO 6000 dock
Device models	Connex VSM 6700 and 6800 series

Ordering Information

06000-900 Braun PRO 6000 dock with thermometer
 106743 Braun PRO 6000 dock (includes rechargeable battery)
 86100 Braun PRO 6000 dock wall-mount kit
 86200 Braun PRO 6000 economy mobile stand-mount kit

Note: Braun ear thermometry is available as a standard option on the Connex Vital Signs Monitor and also as a separately purchased accessory.